



# Dentin bond strength of a universal adhesive system

S. PALMARES<sup>1</sup>, A. PINTO<sup>1</sup>, A. PEQUENO<sup>1</sup>, C. COITO<sup>1</sup>, A. SILVA<sup>1</sup>, R. EIRA. B. SOUSA<sup>1</sup>, A. CAVALHEIRO<sup>1</sup>

<sup>1</sup>University of Lisbon, Lisbon, Portugal

Research performed at UICOB (Biomedical and Oral Sciences Research Unit)

## Purpose

To compare the immediate microtensile bond strength of one universal adhesive system (Scotchbond Universal Adhesive, 3M ESPE, St Paul, MN) applied to dentin according to the etch-and-rinse and the self-etch technique.

## Materials and Methods

Six recently extracted human third molars, intact and without macroscopic evidence of caries or restorations, were assigned to two groups according to the etching strategy: 1) SBU TE D – Scotchbond Universal Adhesive applied as a 2-step etch-and-rinse adhesive on moist dentin and 2) SBU SE D - Scotchbond Universal Adhesive applied as a 1-step self-etch adhesive on moist dentin, both per manufacturer's instructions. Resin composite build-ups (UD4 ENAMEL Plus HRI, Micerium S.p.A. Avegno, GE, Italy) were applied in three increments of 2 mm each, until a height of 6 mm: each layer was light cured for 20 seconds with an additional light polymerization performed on facial, lingual, mesial and distal surfaces for 10 seconds. The teeth were stored in distilled water in an incubator (24h/37°C). Specimens were sectioned to obtain sticks with 1mm<sup>2</sup> of cross sectional area, that were tested to failure in a universal testing machine at a crosshead speed of 1mm/minute, to assess dentin microtensile bond strength ( $\mu$ TBS). A paired-sample t-test was performed when the assumption of normality was valid ( $\alpha=0,05$ ).






Figure 1 - Teeth after being sectioned to obtain sticks.

## Results

The mean  $\mu$ TBS of the SBU SE D group ( $56,9 \pm 2,5$  MPa) was statistically higher than SBU TE D ( $48,0 \pm 2,1$  MPa) ( $p < 0,05$ ).



Graphic 1 - Box-whisker plot of the  $\mu$ TBS for SBU TE D and SBU SE D (x axis represents the group and y axis the MPa).



Graphic 2 - Failure mode distribution: A- adhesive; M – mixed; CD - cohesive in dentin; CC – cohesive in composite.

## Conclusions

It can be concluded that improved bonding effectiveness of Scotchbond Universal Adhesive to dentin seems to be obtained when the adhesive is applied with the self-etch approach.

## References

Hanabusa M, Mine A, Kuboki T, Momoi Y, Van Ende A, Van Meerbeek B, et al. Bonding effectiveness of a new 'multi-mode' adhesive to enamel and dentine. *Journal of dentistry* 2012; 40(6):475-84. Perdigao J, Sezinando A, Monteiro PC. Laboratory bonding ability of a multi-purpose dentin adhesive. *American journal of dentistry* 2012; 25(3):153-8. Mena-Serrano A, Kose C, De Paula EA, Tay LY, Reis A, Loguercio AD, et al. A new universal simplified adhesive: 6-month clinical evaluation. *Journal of esthetic and restorative dentistry : official publication of the American Academy of Esthetic Dentistry* [et al] 2013; 25(1):55-69. Munoz MA, Luque I, Hass V, Reis A, Loguercio AD, Bombarda NH. Immediate bonding properties of universal adhesives to dentine. *Journal of dentistry* 2013; 41(5):404-11. Perdigao J, Kose C, Mena-Serrano A, De Paula E, Tay L, Reis A, et al. A New Universal Simplified Adhesive: 18-Month Clinical Evaluation. *Operative dentistry* 2013a. Marchesi G, Frassetto A, Mazzoni A, Apolonio F, Diolosa M, Cadenaro M, et al. Adhesive performance of a multi-mode adhesive system: 1-Year *in vitro* study. *Journal of dentistry* 2013. Wagner A, Wendler M, Petschelt A, Belli R, Lohbauer U. Bonding performance of universal adhesives in different etching modes. *Journal of dentistry* 2014.